

## Cynomolgus Siglec-15/CD33L3 protein (Recombinant) (C-His) (STJP001292)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, 0.1M Glycine, pH 9.0.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Immunogen Region	Phe20-Thr263
Immunogen Sequence	FVRTKIDTTENLLNTEVHSS PAQRWSMQVPAEVSAAAGDA AVLPCFTFTHPHRHYDGPLTA IWRAGEPYAGPQVFRCAAR GSELCQTALSLHGRFRLGN PRRNDLSLRVERLALADDDR YFCRVEFAGDVH DRYESRHG VRLHVTAAPRIINISVLP GP AHAFRALCTAEGEPPALAW SGPALGN GSAAVPSSGQGHG HLVTAE L PALNHDGRYTCTA ANSLGRSEASVYLF RFHGA

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

