

HCoV-OC43 Hemagglutinin esterase/HE protein (Recombinant) (C-His) (STJP001295)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 µg/ml filtered solution in PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Leu17-Val392
Immunogen Sequence	LGFYNPPTNNVSHVNGDWFL FGDSRSDCNHIVNINPHNYS YMDLNPVLCDSGKISSKAGN SIFRSFHFTDFYNYTGEGQQ IIFYEGVNFTPYHAFKCNRS GSNDIWMQNKGLFYTQVYKN MAVYRSLTFVNVPPYVYNGSA QATALCKSGSLVLNNPAYIA PQANS GDYYYKVEADFYLSG CDEYIVPLCIFNGKFLSNTK YYDDSQYYFNKDTGVYIYGLN STETITTGFDLNCYYLVLP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

