

## Mouse Noggin/NOG protein (Recombinant) (C-His) (STJP001303)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Formulation         | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.                                 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |  |
|--------------------|--|
| Gene ID            | <a href="#">18121</a>  |
| Gene Symbol        | <a href="#">Nog</a>  |
| UniProt ID         | <a href="#">NOGG_MOUSE</a>   |
| Immunogen Region   | Gln28-Cys232   |
| Immunogen Sequence | QHYLHIRPAPSDNPLVLDLI EHPDPIFDPKKDLNETLLR SLLGGHYDPGFMATSPPEDR<br>PGGGGPAGGAEDLAELDQL LRQRPSGAMPSEIKGLEFSE GLAQQKKQRLSKKLRRKLMQ<br>WLWSQTFCPVLYAWNLDLGSR FWPRYVKVGSCFSKRSCSVP EGMVCKPSKSVHLTVLRWRC<br>QRRGGQRCGWIPQIPIISE CKCSC |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

