

Human LILRA5/CD85f protein (Recombinant) (C-His) (STJP001310)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	353514
Gene Symbol	LILRA5
UniProt ID	LIRA5_HUMAN
Immunogen Region	Gly42-Arg268
Immunogen Sequence	GNLSKATLWAEPGSVISRGN SVTIRCQGTLEAQEYRLVKE GSPEPWDTQNPLEPKNKARF SIPSMTEHHAGRYRCYYYS PAGWSEPSDPLELVVTGFYN KPTLSALPSPVVTSGENVTL Q CGSRLRFDRFILTEEGDH KLSWTLDSQLTPSGQFQALF PVGPVTPSHRWMLRCYGSRR HI LQVWSEPSDLLEIPVSG AADNLSPSQNKSDSGTASHL QDYAVENLIR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

