

## Human FLT-3/FLK-2 protein (Recombinant) (C-hFC&His) (STJP001311)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">2322</a>
Gene Symbol	<a href="#">FLT3</a>
UniProt ID	<a href="#">FLT3_HUMAN</a>
Immunogen Region	Asn27-Asn541
Immunogen Sequence	NQDLPIKCVLINHKNDSS V GKSSSYPMVSESPEDLGCA LRPQSSGTVYEAAAVEVDVS ASITLQVLVDAPGNISCLWV FKHSNLNCQPHFDLQNRGVV SMVILKMTETQAGEYLLFIQ SEATNYTILFTVSIRNTLLY TLRRPYFRKMENQDALVCIS ESVPPIVEWVLCDSQGESC KEESPAVVKKEEKVLHELFG TDIRCCARNELGRECTRLFT IDLNQTPQTTLPQLFLKVG

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

