

Rat VEGF-A/VEGF164 protein (Recombinant) (C-His) (STJP001318)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	83785
Gene Symbol	Vegfa
UniProt ID	VEGFA_RAT
Immunogen Region	Ala27-Arg190
Immunogen Sequence	APTTEGEQKAHEVVKFMDVY QRSYCRPIETLVDFQEYPD EIEYIFKPSCVPLMRCAGCC NDEALECVPTSESNTMQIM RIKPHQSQHIGEMSFLLQHSR CECRPKKDRTPENHCEPCS ERRKHLFVQDPQTCKCCKN TDSRCKARQLELNERTCRCD KPRR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

