

Mouse CTLA-4/CD152 protein (Recombinant) (C-His) (STJP001319)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12477
Gene Symbol	Ctla4
UniProt ID	CTLA4_MOUSE
Immunogen Region	Glu36-Asp161
Immunogen Sequence	EAIQVTQPSVVLASSHGVS FPCEYSPSHNTDEV RVTVLR QTNDQMTEVCATTFTTEKNTV GFLDYPPFCSGTFNFSRVNL TIQGLRAVD TGLYLCKVELM YPPPYFVG MNGTQIYVIDP E PCPDSD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

