

Human LILRA3/CD85e protein (Recombinant) (C-His) (STJP001323)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11026
Gene Symbol	LILRA3
UniProt ID	LIRA3_HUMAN
Immunogen Region	Gly24-Glu439
Immunogen Sequence	GPLPKPTLWAEPGSGVITQGS PVTLCQGSLETQEYHLYRE KKTALWITRIPQELVKKGQF PILSITWEHAGRYCCYGS HTAGLSESSDPLELVVTGAY SKPTLSALPSPVVTSGGNVT I QCDSQVAFDGFILCKEKE DEHPQCLNSHSHARGSSRAI FSVGVPVSPRRWSYRCYGYD SR APYVWSLPSDLLGLLVP GVSKKPSLSVQPGPVVAPGE KLTFQCGSDAGYDRFVLYKE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

