

Human TIMP-1 protein (Recombinant) (C-His) (STJP001333)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7076
Gene Symbol	TIMP1
UniProt ID	TIMP1_HUMAN
Immunogen Region	Cys24-Ala207
Immunogen Sequence	CTCVPPHPQTAFCSNDLVIR AKFVGTPENVNQTTLQRYEI KMTKMYKGFQALGDAADIRF VYTPAMESVCGYFHRSHNRS EEFLIAGKLQDGLLHITCS FVAPWNSLSLAQRRGFTKTY TVGCEECTVFPCLSLPCKLQ SGTHCLWTDQLLQGSSEKGFQ SRHLACLPREPLCTWQSLR SQIA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

