

Human Prolyl endopeptidase FAP protein (Active Recombinant) (C-His) (STJP001339)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Supplied as a 0.22 Mu m filtered solution in PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2191
Gene Symbol	FAP
UniProt ID	SEPR_HUMAN
Immunogen Region	Leu26-Asp760
Immunogen Sequence	LRPSRVHNSEENTMRALTLK DILNGTFSYKTFPPNWISGQ EYLHQSADNNIVLYNIETGQ SYTILSNRTMKSVNASNYGL SPDRQFVYLES DYSKLWRY S YTATYYIYDLSNGEFVRGNE LPRPIQYLCWSPVGSKLAYV YQNNIYKQRPDPPFQITF NGRENKIFNGIPDWVVEEEM LATKYALWWSPNGKFLAYAE FNDTIPVIAYSYYGDEQYP RTINIPYPKAGAKNPVVR I

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

