

## Human Dkk-1 protein (Recombinant) (C-His) (STJP001340)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">22943</a>
Gene Symbol	<a href="#">DKK1</a>
UniProt ID	<a href="#">DKK1_HUMAN</a>
Immunogen Region	Thr32-His266
Immunogen Sequence	TLNSVLNSNAIKNLPPLGG AAGHPGSAVSAAPGILYPGG NKYQTDINYQPYPCAEDEEC GTDEYCASPTRGGDAGVQIC LACRKRKRCMRHAMCCPGN YCKNGICVSSDQNHFRGEIE ETITESFGNDHSTLDGYSRR TTLSSKMYHTKGQEGSVCLR SSDCASGLCCARHFWSKICK PVLKEGQVCTKRRRKGSHGL EIFQRCCYCGEGLSRIQKDH HQASNSSLRLHTCQRH

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

