

Human SIRP-beta-1/SIRPB1/CD172b protein (Recombinant) (C-His) (STJP001343)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10326
Gene Symbol	SIRPB1
UniProt ID	SIRB1_HUMAN
Immunogen Region	Glu30-Leu371
Immunogen Sequence	EDELQVIQPEKSVSVAAGES ATLRCAMTSLIPVGPIMWFR GAGAGRELIYNQKEGHFPRV TTVSELTKRNNLDFSISIN ITPADAGTYCVKFRKGGSPD DVEFKSGAGTELSVRAKPSA PVVSGPAVRATPEHTVSFTC ESHGFSPRDITLKWFKNGNE LSDFQTNVDPAGDSVSYSIH STARVVLTRGDVHSQVCEI AHITLQGDPLRGTANLSEAI RVPPTLEVTQQPMRAENQA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

