

Human TNFRSF21/DR6 protein (Recombinant) (C-hFc&His) (STJP001345)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27242
Gene Symbol	TNFRSF21
UniProt ID	TNR21_HUMAN
Immunogen Region	Gln42-Leu350
Immunogen Sequence	QPEQKASNLIGTYRHVDRAT GQVLTCDKCPAGTYVSEHCT NTSLRVCSSCPVGTFTTRHEN GIEKCHDCSQPCPWPMIEKL PCAALTDRECTCPPGMFQSN ATCAPHTVCPVGGWVRRKKG ETEDVRCKQCARGTFSDVPS SVMKCKAYTDCLSQLLVVIK PGTKETDNVCGTLPFSST SPSPGTAIFRPEHMETHEV PSTYVPKGMNSTESNSSAS VRPKVLSSIQEGTVPDNTS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

