

## Human CD34 protein (Recombinant) (C-hFc&His) (STJP001349)

### GENERAL INFORMATION

|               |              |
|---------------|--------------|
| Product Type  | Proteins     |
| Host / Source | HEK293 cells |

### PRODUCT PROPERTIES

|                     |                                                                                                |
|---------------------|------------------------------------------------------------------------------------------------|
| Formulation         | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.                                 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                                  |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">947</a>                                                                                                                                                                                                                                              |
| Gene Symbol        | <a href="#">CD34</a>                                                                                                                                                                                                                                             |
| UniProt ID         | <a href="#">CD34_HUMAN</a>                                                                                                                                                                                                                                       |
| Immunogen Region   | Ser32-Thr290                                                                                                                                                                                                                                                     |
| Immunogen Sequence | SLDNNGTATPELPTQGTFN VSTNVSYQETTPSTLGSTS LHPVSQHGNEATTNITETTV<br>KFTSTSVITSVYGNTNSSVQ SQTSVISTVFTTPANVSTPE TTLKPSPGPNVSDLTSTTST<br>SLATSPTKPYTSSSPILSDI KAEIKCSGIREVKLTQGICL EQNKTSSCAEFKKDRGEGLA<br>RVLCGEEQADADAGAVCSL LLAQSEVRPQCLLLVLANRT EISSKLQLMKKHQSDLKKL |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

