

Human B7-H6/NCR3LG1 protein (Recombinant) (C-His) (STJP001350)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	374383
Gene Symbol	NCR3LG1
UniProt ID	NR3L1_HUMAN
Immunogen Region	Asp25-Ser262
Immunogen Sequence	DLKVEMMAGGTQITPLNDNV TIFCNIFYSQPLNITSMGIT WFWKSLTFDKEVKVFEFFGD HQEAFRPGAIVSPWRLKSGD ASLRLPGIQLEEEAGEYRCEV VVTPLKAQGTVQLEVVASPA SRLLDQVGMKENEDKYMCE SSGFYPEAINITWEKQTQKF PHPIEISEDVITGPTIKNMD GTFNVTSCCLKLNSSQEDPGT VYQCVVRHASLHTPLRSNFT LTAARHSLSETEKTDNFS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

