

Human IL-3 protein (Active Recombinant) (C-His) (STJP001354)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3562
Gene Symbol	IL3
UniProt ID	IL3_HUMAN
Immunogen Region	Ala20-Phe152
Immunogen Sequence	APMTQTTPKTSWVNCNMI DEIITHLKQPPLPLDFNNL NGEDQDILMENNLRPPNLEA FNRAVKSLQNASAIESILKN LLPCLPLATAAPTRHPIHK DGDWNEFRKLTFLYKLTLEN AQAQQTLSLAIF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

