

Human IGFBP-3 protein (Active Recombinant) (C-His) (STJP001355)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3486
Gene Symbol	IGFBP3
UniProt ID	IBP3_HUMAN
Immunogen Region	Gly28-Lys291
Immunogen Sequence	GASSAGLGPVVRCEPCDARA LAQCAPPPAVCAELVREPGC GCCLTCALSEGQPCGIYTER CGSGLRCQSPDEARPLQAL LDGRGLCVNASAVSRLRAYL LPAPPAPGNASESEEDRSAG SVESPSVSSTHRVSDPKFHP LHSKIIIIKKGHAKDSQRYK VDYEQSTDTQNFSSSESKRE TEYGPCRREMEDTLNHLKFL NVLSPRGVHIPNCDKKGFYK KKQCRPSKGRKRGFCWCVD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

