

Human Complement factor D/CFD protein (Recombinant) (C-His) (STJP001360)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1675
Gene Symbol	CFD
UniProt ID	CFAD_HUMAN
Immunogen Region	Ile26-Ala253
Immunogen Sequence	ILGGREAEAHARPYMASVQL NGAHLCCGGVLVAEQWVLSAA HCLEDAADGKQVLLGAHSL SQPEPSKRLYDVLRAVPHPD SQPDTIDHLLLLQLSEKAT LGPAVRPLPWQRVDRDVAPG TLCDVAGWGI VNHAGRRPDS LQHVLLPVLDRA TCNRRRTHH DGAITERLMCAESNRRD SCK GDSGGPLVCGGVLEGVVTSG SRVCGNRKKPGIYTRVASYA AWIDSVLA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

