

Human TIMP-2 protein (Recombinant) (C-His) (STJP001362)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7077
Gene Symbol	TIMP2
UniProt ID	TIMP2_HUMAN
Immunogen Region	Cys27-Pro220
Immunogen Sequence	CSCSPVHPQQAFCNADVIR AKAVSEKEVDSGNDIYGNI KRIQYEIKQIKMFKGPEKDI EFIYTAPSSAVCGVSLDVGG KKEYLIAGKAEGDGKMHITL CDFIVPWTDLSTTQKKSLNH RYQMGCECKITRCMPICYI SSPDECLWMDWVTEKNINGH QAKFFACIKRSDGSCAWYRG AAPPKQEFLDIEDP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

