

## Human LILRB2/ILT-4/CD85d protein (Recombinant) (C-hFc&His) (STJP001364)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">10288</a>
Gene Symbol	<a href="#">LILRB2</a>
UniProt ID	<a href="#">LILRB2_HUMAN</a>
Immunogen Region	Gln22-Val461
Immunogen Sequence	QTGTIPKPTLWAEPSVITQ GSPVTLSCQGSLEAQEYRLY REKKSASWITRIRPELVKNG QFHIPSITWEHTGRYGCQYY SRARWSELSDPLVLMVTGAY PKPTLSAQSPVVTSGGRVT LQCESQVAFGGFILCKEGEE EHPQCLNSQPHARGSSRAIF SVGPVSPNRRWHRWCYGYDL NSPYVWSSPSDLELLVPGV SKKPSLSVQPGPVPVAPGESL TLQCVSDVGYDRFVLYKEG

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

