

Mouse IL-15RA/CD215 protein (Recombinant) (C-hFc&His) (STJP001365)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16169
Gene Symbol	Il15ra
UniProt ID	I15RA_MOUSE
Immunogen Region	Gly33-Lys205
Immunogen Sequence	GTTCPPPVSIHADIRVKNY SVNSRERYVCNSGFKRKAGT STLIECVINKNTNVAHWTPP SLKCIRDPSLAHYSPVPTVV TPKVTSQPESPPSAKEPEA FSPKSDTAMTTETAIMPGSR LTPSQTTAGTTGTGSHKSS RAPSLAATMTLEPTASTSLR ITEISPHSSKMTK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

