

Human Frizzled-4/FZD4/CD344 protein (Recombinant) (C-hFc&His) (STJP001367)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8322
Gene Symbol	FZD4
UniProt ID	FZD4_HUMAN
Immunogen Region	Phe37-Glu180
Immunogen Sequence	FGDEEERRCDPIRISMQNL GYNVTKMPNLVGHELQTDAE LQLTTFTPLIQGCSSQLQF FLCSVYVPMCTEKINIPGP CGGMCLSVKRRCEPVLKEFG FAWPESLNCSKFFPPQNDHNH MCMEGPGDEEVPLPHKTIQ PGEE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

