

Human VEGFR-2/KDR/CD309 protein (Recombinant) (C-hFc&His) (STJP001370)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3791
Gene Symbol	KDR
UniProt ID	VGFR2_HUMAN
Immunogen Region	Ala20-Glu764
Immunogen Sequence	ASVGLPSVSLDLPRLSIQKD ILTIKANTTLQITCRGQDDL DWLWPNNQSGSEQRVEVTEC SDGLFCKTLTIPKVIGNDTG AYKCFYRETDLASVIYVYVQ DYRSPFIASVSDQHGVVYIT ENKNKTVVIPCLGSISNLNV SLCARYPEKRFVDPGNRIW DSKKGFTIPSYMISYAGMVF CEAKINDESYQSIMYIVVVV GYRIYDVVLSPSHGIELSVG EKLVLNCTARTELVGIDF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

