

Human Vitronectin/V75/VTN protein (Recombinant) (C-His) (STJP001377)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7448
Gene Symbol	VTN
UniProt ID	VTNC_HUMAN
Immunogen Region	Asp20-Leu478
Immunogen Sequence	DQESCKGRCTEGFNVDKCKQ CDELCSYYQSCCTDYTAECK PQVTRGDVFTMPED EYTVYD DGEEKNNATVHEQVGGPSLT SDLQAQSKGNPEQTPVLKPE EEAPAPEVGASKPEGIDSRP ETLHPGRPQPPAEEELCSGK PFDAFTDLKNGSLFAFRGQY CYELDEKAVRPGYPKLRDV WGIEGPIDAAFTRINCQGKT YLFKGSQYWR FEDGVLD PDY PRNISDGFDPDNVDAAL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

