

Human TNFRSF10B/DR5/TRAIL-R2 protein (Recombinant) (C-hFc&His) (STJP001381)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20Â°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8795
Gene Symbol	TNFRSF10B
UniProt ID	TR10B_HUMAN
Immunogen Region	Ile56-Glu182
Immunogen Sequence	ITQQDLAPQQRAAPQQKRSS PSEGLCPPGHHISEDGRDCI SCKYGQDYSTHWNDLLFCLR CTRCDSGEVELSPCTTTRNT VCQCEEGTFREEDSPEMCRK CRTGCPRGMVKVGDCTPWSD IECVHKE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

