

Mouse Fc-gamma RIV/CD16-2 protein (Recombinant) (C-His) (STJP001383)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	246256
Gene Symbol	Fcgr4
UniProt ID	FCG3A_MOUSE
Immunogen Region	Gly21-Gln 203
Immunogen Sequence	GLQKAVVNLDPKWVRVLEED SVTLRCQGTFSPEDNSIKWF HNESLIPHQDANYVIQSARV KDSGMYRCQTALSTISDPVQ LEVHMGWLLLQTTKWLFEQEG DPIHLRCHSWQNRPVVKVTY LQNGKGKKYFHENSELLIPK ATHNDSGSYFCRGLIGHNNK SSASFRISLGDPGSPSMFPP WHQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

