

Human TNFRSF18/CD357 protein (Recombinant) (C-hFc&His) (STJP001388)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8784
Gene Symbol	TNFRSF18
UniProt ID	TNR18_HUMAN
Immunogen Region	Gln26-Glu161
Immunogen Sequence	QRPTGGPGCGPGRLLLTGT DARCCRVHTTRCCRDYPGEE CCSEWDCMCVQPEFHCGDPC CTTCRHHPGPPGQGVQSQK FSGFQCIDCASGTFSGGHE GHCKPWTDCTQFGFLTVPFG NKTHNAVCPGSPPAE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

