

Human CCL21 protein (Recombinant) (No tag) (STJP001396)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6366
Gene Symbol	CCL21
UniProt ID	CCL21_HUMAN
Immunogen Region	Ser24-Pro134
Immunogen Sequence	SDGGAQDCCLKYSQRKIPAK VVRSYRKQEPSLGCIPAIL FLPRKRSQAELCADPKELWV QQLMQHLDKTPSPQKPAQGC RKDRGASKTGKKGKGSKGCK RTERSQTPKGP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

