

Human CDH6/K-Cadherin protein (Recombinant) (C-hFc&His) (STJP001398)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1004
Gene Symbol	CDH6
UniProt ID	CADH6_HUMAN
Immunogen Region	Thr19-Ala615
Immunogen Sequence	TLSTPLSKRTSGFPAKKRAL ELSGNSKNELNRSKRSMWN QFFLLEEYTGSDYQYVGKLH SDQDRGDGSLKYILSGDGAG DLFIIINENTGDIQATKRLDR EEKPVYILRAQAINRRTGRP VEPESEFIIKIHINDNEPI FTKEVYTATVPMSDVGTFV VQVTATDADDPTYGNSAKVV YSILQGQPYFSVESETGIK TALLNMDRENREYQVVIQA KDMGGQMGGLSGTTTVNIT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

