

Human Fc epsilon RII/CD23 protein (Recombinant) (N-his) (STJP001400)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2208
Gene Symbol	FCER2
UniProt ID	FCER2_HUMAN
Immunogen Region	Asp 48-Ser321
Immunogen Sequence	DTTQSLKQLEERAARNVSQV SKNLESHHGDMQAQKSQSTQ ISQELEELRAEQRLKSQDL ELSWNLNGLQADLSSFKSQE LNERNEASDLLERLREEVTK LRMELVSSGFVNCNTCPEKW INFQRKCYFFGKGTKQWVHA RYACDDMEGQLVSIHSPEEQ DFLTKHASHTGSWIGLRNLD LKGFIWVDGSHVDYSNWAP GEPTSRSQGEDCVMMRGSGR WNDAFCDRKLGAWVCDRLA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

