

## Human CCL2/MCP-1 protein (Recombinant) (No tag) (STJP001403)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6347</a>
Gene Symbol	<a href="#">CCL2</a>
UniProt ID	<a href="#">CCL2_HUMAN</a>
Immunogen Region	Gln24-Thr99
Immunogen Sequence	QPDAINAPVTCCYNFTNRKI SVQRLASYRRITSSKCPKEA VIFKTIVAKEICADPKQKWW QDSMDHLDKQTQTPKT

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

