

Human Endoglin/CD105 protein (Recombinant) (C-His) (STJP001405)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2022
Gene Symbol	ENG
UniProt ID	EGLN_HUMAN
Immunogen Region	Glu26-Gly586
Immunogen Sequence	ETVHCDLQPVGPERGEVTTYT TSQVSKGCVAPNAILEVH VLFLEFPTGPSQLELTQAS KQNGTWPREVLLVLSVNSSV FLHLQALGIPLHLAYNSSLV TFQEPPGVNTELPSPFKTQ ILEWAAERGPITSAAELNDP QSILLRLGQAQGSLSFCMLE ASQDMGRTLWRPRTPALVR GCHLEGVAGHKEAHILRVLP GHSAGPRTVTVKVELSCAPG DLDVAVLILQGPPYVSWLID

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

