

Human VEGF-D/FIGF protein (Recombinant) (C-His) (STJP001416)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 M filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2277
Gene Symbol	VEGFD
UniProt ID	VEGFD_HUMAN
Immunogen Region	Phe93-Ser201
Immunogen Sequence	FYDIETLKVIDEEWQRTQCS PRETCVEVASELGKSTNTFF KPPCVNVFRCGGCCNEESLI CMNTSTSYISKQLFEISVPL TSVPELVVKVANHTGCKCL PTAPRHPYS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

