

Cynomolgus IL-2RG/CD132 protein (Recombinant) (C-hFc&Avi) (STJP001418)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	Leu23-Asn254
Immunogen Sequence	LNTTILTPNGNEDATTDFFL TSMPTDSLVSSTLPLPEVQC FVFNVEYMNCTWNSSEPQP TNLTLYWYKNSDNDKVQKC SHYLFSEEITSGCQLQEKEI HLYQTFVVQLQDPREPRRQA TQMLKLQNLVIPWAPENLTL RKLSESOLELNWNNRFLNHC LEHLVQYRTDWDHSWTEQSV DYRHKFSLPSVDGQKRYTFR VRSRFNPLCGSAQHWSEWSH PIHWGSNSSKEN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

