

Human Carbonic anhydrase 9 protein (Active Recombinant) (C-His&Avi) (STJP001420)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	768
Gene Symbol	CA9
UniProt ID	CAH9_HUMAN
Immunogen Region	Gln38-Asp414
Immunogen Sequence	QRLPRMQEDSPLGGGSSGED DPLGEEDLPSEEDSPREEDP PGEEDLPGEEDLPGEEDLPE VKPKSEEEGSLKLEDLPTVE APGDPQEPQNNNAHRDKEGDD QSHWRYGGDPPWPRVSPACA GRFQSPVDIRPQLAFCAL RPLELLGFQLPPLPELRLRN NGHSVQLTLPPGLEMALGPG REYRALQLHLHWGAAGRPGS EHTVEGHRFP AEIHVVHLST AFARVDEALGRPGGLAVLA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

