

Human NKAT-1/KIR2DL1/CD158a protein (Recombinant) (C-His&Avi) (STJP001423)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3802
Gene Symbol	KIR2DL1
UniProt ID	K12L1_HUMAN
Immunogen Region	His22-Arg242
Immunogen Sequence	HEGVHRKPSLLAHPGLVKS EETVILQCWSDVMFEHLLH REGMFNDTLRLIGEHHDGVS KANFSISRMTQDLAGTYRCY GSVTHSPYQVSAPSDPLDIV IIGLYEKPSLSAQPGPTVLA GENVTLCSSRSSYDMYHLS REGEAHERRLPAGPKVNGTF QADFPLGPATHGGTYRCFGS FHDSPYEWSKSSDLLSVST GNPSNSWSPSTEPSKGTGNP R

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

