

Human IFN-alpha 2 protein (Recombinant) (C-His) (STJP001424)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3440
Gene Symbol	IFNA2
UniProt ID	IFNA2_HUMAN
Immunogen Region	Cys24-Glu188
Immunogen Sequence	CDLPQTHSLGSRRTLMLLAQ MRKISLFSCLKDRHDFGFPQ EEFGNQFQKAETIPVLHEMI QQIFNLFSTKDSSAAWDETL LDKFYTELYQQLNDLEACVI QGVGVTTETPLMKEDSILAVR KYFQRITLYLKEKKYSPCAW EVVRAEIMRSFSLSTNLQES LRSKE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

