

Human LILRB3/ILT-5/CD85a protein (Recombinant) (C-His&Avi) (STJP001428)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene Symbol	LILRB3 Uniprot ID= LIRB3_HUMAN "
Immunogen Region	Gly24-Glu443
Immunogen Sequence	<p> GPFKPTLWAEPGSVISWGS PVTIWCQGSQEAQEYRLHKE GSPEPLDRNNPLEPKNKARF SIPSMTEHHAGRYRCHYYSS AGWSEPSDPLEMVMTGAYSK PTLALPSPVVASGNNMTRLR CGSQKGYHHFVLMKEGEHQL PRTLDSQQLHSRGFQALFPV GPVTPSHRWRFTCYYYTNT PWVWSPSDPLEILPSGVSR KPSLLTLQGPVLAPGQSLTL QCGSDVGYNRFVLYKEGER </p>

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

