

Human FGF-6 protein (Recombinant) (No tag) (STJP001429)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	E.coli

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2251
Gene Symbol	FGF6
UniProt ID	FGF6_HUMAN
Immunogen Region	Gly41-Ile208
Immunogen Sequence	GTRANNTLLDSRGWGTLLSR SRAGLAGEIAGVNWESGYLV GIKRQRRLYCNVGIGFHLQV LPDGRISGTHEENPYSLEI STVERGVVSLFGVRSALFVA MNSKGRLYATPSFQEECKFR ETLLPNNYNAYESDLYQGTY IALSKYGRVKRGSKVSPIMT VTHFLPRI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

