

Human LAIR-1/CD305 protein (Recombinant) (C-His&Avi) (STJP001432)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3903
Gene Symbol	LAIR1
UniProt ID	LAIR1_HUMAN
Immunogen Region	Gln22-His163
Immunogen Sequence	QEEDLPRPSISAEPGTVIPL GSHVTFVCRGPVGVQTFRLE RESRSTYNDTEDVSQASPSE SEARFRIDSVSEGNAGPYRC IYYKPPKWSEQSDYLELLVK ETSGGPDSPDTEPGSSAGPT QRPSDNSHNEHAPASQGLKA EH

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

