

Human Apolipoprotein A-II/APOA2 protein (Recombinant) (C-hFc&Avi) (STJP001433)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	336
Gene Symbol	APOA2
UniProt ID	APOA2_HUMAN
Immunogen Region	Ala19-Gln100
Immunogen Sequence	ALVRRQAKEPCVESLVSQYF QTVTDYDGKDLMEKVKSPQLQ AEAKSYFEKSKEQLTPLIKK AGTELVNFLSYFVELGTQPA TQ

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

