

Human Apolipoprotein E/ApoE protein (Recombinant) (C-His) (STJP001435)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	348
Gene Symbol	APOE
UniProt ID	APOE_HUMAN
Immunogen Region	Lys19-His317
Immunogen Sequence	KVEQAVETEPEPELRQQTEW QSGQRWELALGRFWDYLRWV QTLSEVQVEELLSSQVTQEL RALMDETMKELKAYKSELEE QLTPVAEETRALRSKELQAA QARLGADMEDVCGRLVQYRG EVQAMLGQSTEELRVRLASH LRKLRKRLLRDADDLQKRLA VYQAGAREGAERGLSAIRER LGPLVEQGRVRAATVGSAG QPLQERAQAWGERLRARMEE MGSRTDRDLDEVKEQVAEV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

