

Human gp140/CDCP1/CD318 protein (Recombinant) (C-His) (STJP001442)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64866
Gene Symbol	CDCP1
UniProt ID	CDCP1_HUMAN
Immunogen Region	Phe30-Leu666
Immunogen Sequence	FEIALPRESNITVLIKLGTP TLLAKPCYIVISKRHITMLS IKSGERIVFTFSCQSPENHF VIEIQKNIDCMGSPCFGEV QLQPSTSLPTLNRTFIWDV KAHKSIGLELQFSIPRLRQI GPGESCPDGVTHSISGRIDA TVVRIGTFCSNGTVSRIKMQ EGVKMAHLFPWFHPRNVSGF SIANRSSIKRLCIIESVFEG EGSATLMSANYPEGFPEDEL MTWQFVVPAPHLRASVSFLN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

