

Human Mature TGF-beta 1 protein (Active Recombinant) (No tag) (STJP001449)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of 50 mM Glycine, 150 mM NaCl, pH 3.5.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7040
Gene Symbol	TGFB1
UniProt ID	TGFB1_HUMAN
Immunogen Region	Ala279-Ser390
Immunogen Sequence	LSTSKTIDMELVKRKRIEAI RGQLSKLRLASPPSQGEVP PGPLPEAVLALYNSTRDRVA GESAEPEPEPEADYYAKEVT RVLMMVETHNEIYDKFKQSTH SIYMFNTSELREAVPEPVL LSRAELRLLRLKLVKVEQHVE LYQKYSNNSWRYLSNRLAP SDSPEWLSFDVTGVVQWLS RGGEIEGFRLSAHCSDSRD NTLQVDINGFTTGRRGLAT IHGMNRPFLMATPLERA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

