

Human Adiponectin/Acrp30/ADIPOQ protein (Recombinant) (C-His) (STJP001454)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9370
Gene Symbol	ADIPOQ
UniProt ID	ADIPO_HUMAN
Immunogen Region	Glu19-Asn244
Immunogen Sequence	ETTTQGGVLLPLPKGACTG WMAGIPGHPGHNGAPGRDGR DGTPGEKGEKGDPLIGPKG DIGETGVPGAEGPRGFPGIQ GRKGEPGEGAYVYRSAFEVSG LETYVTIPNMPIRFTKIFYN QQNHYDGGSTGKFCNIPGLY YFAYHITVYMKDVKVSFLFKK DKAMLFYDQYQENNVQAS GSVLLHLEVGQWVLQVYGE GERNGLYADNDNDSTFTGFL LYHDTN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

