

Mouse Adiponectin/Acrp30/ADIPOQ protein (Recombinant) (C-His) (STJP001459)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11450
Gene Symbol	Adipoq
UniProt ID	ADIPO_MOUSE
Immunogen Region	Glu18-Asn247
Immunogen Sequence	EDDVTTEELAPALVPPPKG TCAGWMAGIPGHPGHNGTPG RDGRDGTGPGEKGEKGDAGLL GPKGETGDVGMTGAEGPRGF PGTPGRKGEPGEAAVYVYRSA FSVGLETRVTVPNVPIRFTK IFYNQNNHYDGGSTGKFYCN I PGLYYFSYHITVYMKDVKVS LFKKDKAVLFYDQYQEKNV DQASGSVLLHLEVGDVWLQ VYGDGDHNGLYADNVNDSTF TGFLLYHDTN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

