

Mouse IL-2-RA/CD25 protein (Recombinant) (C-hFc) (STJP001474)

GENERAL INFORMATION

| | |
|---------------|--------------|
| Product Type | Proteins |
| Host / Source | HEK293 cells |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Formulation | Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4. |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|--------------------|---|
| Gene ID | 16184 |
| Gene Symbol | Il2ra |
| UniProt ID | IL2RA_MOUSE |
| Immunogen Region | Glu22-Lys236 |
| Immunogen Sequence | ELCLYDPPEVPNATFKALSY KNGTILNCECKRGFRRLKEL VYMRCLGNSWSSNCQCTSNS HDKSRKQVTAQLEHQKEQQT TTDMMQKPTQSMHQENLTGHC REPPPWKHEDSKRIYHFVEG QSVHYECIPGYKALQRGPAI SICKMKCGKTGWTQPQLTCV DEREHHRFLASEESQGSFRNS SPESETSCPITTTDFPQTE TTAMTETFVLTMEYK |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

