

## Human CSF1R/M-CSF R/CD115 protein (Recombinant) (C-hFc&His) (STJP001477)

### GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

### PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1436</a>
Gene Symbol	<a href="#">CSF1R</a>
UniProt ID	<a href="#">CSF1R_HUMAN</a>
Immunogen Region	Ile20-Glu512
Immunogen Sequence	IPVIEPSVPELVVKGATVT LRCVNGSVVEWDGPPSPHWT LYSDGSSSILSTNNATFQNT GTYRCTEPGDPLGGSAIHL YVKDPAWPWNVLAQEVVFE DQDALLPCLLTDPVLEAGVS LVRVGRPLMRHTNYSFSPW HGFTIHRAKFIQSQDYQCSA LMGGRKVMISIRLKVQKVI PGPPALTLPVPAELVRIRGEA AQIVCSASSVDVNFDFLQH NNTKLAIPQQSDFHNNRYQ

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

