

Mouse Siglec-3/CD33 protein (Recombinant) (C-His) (STJP001494)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12489
Gene Symbol	Cd33
UniProt ID	CD33_MOUSE
Immunogen Region	Asp18-Glu240
Immunogen Sequence	DLEFQLVAPESVTVEEGLCV HVPCSVFYPSIKLTGCVTG SWLRKGVSLHEDSPVATSDP RQLVQKATQGRFQLLGDPQK HDCSLFIRDAQKNDTGMYFF RVVREPFRVRSYKKSQSLSLH VTLSRTPDIIIPGTLEAGY PSNLTCSVPWACEQGTPTTF SWMSTALTSLSRRTDSSVL TFTPQPQDHGKLTCLVTFS GAGVTVERTIQLNVTRKSGQ MRE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

