

Mouse c-Kit/CD117 protein (Recombinant) (C-His) (STJP001497)

GENERAL INFORMATION

Product Type	Proteins
Host / Source	HEK293 cells

PRODUCT PROPERTIES

Formulation	Lyophilised from a 0.22 Mu m filtered solution of PBS, pH 7.4.
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16590
Gene Symbol	Kit
UniProt ID	KIT_MOUSE
Immunogen Region	Ser25-Pro527
Immunogen Sequence	SQPSASPGEPSPPSIHPAQS ELIVEAGDTLSLTCIDPDFV RWTFKTYFNEMVENKKNEWI QEKAETRGTGTYTCSNSNGL TSSIYVVRDPAKLFLVGLP LFGKEDSDALVRCPLTDPQV SNYSLECDGKSLPTDLTFV PNPKAGITIKNVKRAYHRLC VRCAAQRDGTWLHSDKFTLK VRAAIKAI PVVSV PETSHLL KKGDTFTVVCTIKDVSTSVN SMWLKMN PQQHIAQVKHN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

